



# GSS23GSK/GMK

## GE Appliances 23.2 Cu. Ft. Side-By-Side Refrigerator

### DIMENSIONS AND INSTALLATION INFORMATION (IN INCHES)

OVERALL DIMENSIONS	Height to top of hinge (in.) <b>A</b>	69-1/2
	Height to top of case (in.) <b>B</b>	69
	Case depth without door (in.) <b>C</b>	28-1/2
	Case depth less door handle (in.) <b>D</b>	32-3/4
	Case depth with door handle (in.) <b>E</b>	35-1/4
	Depth with fresh food door open 90° (in.) <b>F</b>	50-1/2
	Width (in.) <b>G</b>	32-3/4
	Width with door open 90° with door handle (in.) <b>H</b>	38-1/4
AIR CLEARANCES	Each side (in.)	1/8
	Top (in.)	1
	Back (in.)	2



If installed against a wall, allow clearance of 14-3/16" on freezer side to remove bin.

Clearance required to remove fresh food full-size pan without disassembling is 18-15/16".

Allow additional space for any necessary leveling adjustments.



For answers to your Monogram, GE Café™ Series, GE Profile™ Series or GE Appliances product questions, visit our website at [geappliances.com](http://geappliances.com) or call GE Answer Center® Service, 800.626.2000.



Specification Revised 9/16

# GSS23GSK/GMK

## GE Appliances 23.2 Cu. Ft. Side-By-Side Refrigerator

### FEATURES AND BENEFITS

Integrated shelf support system – Provides strong, flexible support

Adjustable slide-out, spillproof glass shelves – Raised edges help contain spills and make cleanup quick and easy

Fresh food multi-level drawers – Provide the ideal environment for storing your fruits and vegetables

Adjustable gallon door bins – Offer ideal space for storing gallon containers in the door

Arctica icemaker – A special access door puts ice within reach

Advanced water filtration uses MWF replacement filter – Reduces trace pharmaceuticals from water and ice\* (\* Removes 98% of ibuprofen, atenolol, fluoxetine, progesterone and trimethoprim. These pharmaceuticals are not necessarily in all users' water)

Model GSS23GSKSS – Stainless steel

Model GSS23GMKES – Slate

